


BREATHABLE BACK UPHOLSTERY

TENSION ADJUSTABLE FOR TiLite™ ZR (SERIES 2), ZRA (SERIES 2), AERO Z (SERIES 2), ZR, ZRA, AERO Z, Z10, TR, YR, BB, AERO R, TR (SERIES 3), AERO T, TRA (SERIES 2)



ITEM	QUANTITY	PART NO.	DESCRIPTION
1	1	TABLE THREE	BACKREST EAR RIGHT (SEE NEXT PAGE)
2	1	TABLE THREE	BACKREST EAR LEFT (SEE NEXT PAGE)
3	1	TABLE FOUR	BACK PAD TENSION ADJUSTABLE (SEE NEXT PAGE)
4*	1-3	TABLE TWO	BACK STRAP (SEE NEXT PAGE)
1-4**		TABLE ONE	TENSION ADJUSTABLE, BACK UPHOL. ASSEMBLY W/ PUSH (SEE BELOW)
1-4**		TABLE ONE	TENSION ADJUSTABLE, BACK UPHOLSTERY ASSEMBLY (SEE BELOW)

TABLE ONE - TENSION ADJUSTABLE BACKREST ASSEMBLIES

CHAIR WIDTH	STANDARD ADJ. HT 8.5"-11.5"	STANDARD ADJ. HT 12"-16.5"	STANDARD ADJ. HT 16"-20.5"	ADJ. HT 10"-12.5" W/PUSH HANDLES	ADJ. HT 12"-16.5" W/PUSH HANDLES	ADJ. HT 16"-20.5" W/PUSH HANDLES
10"	AS115332	AS115343	AS115354	AS115365	AS115376	AS115387
11"	AS115333	AS115344	AS115355	AS115366	AS115377	AS115388
12"	AS115334	AS115345	AS115356	AS115367	AS115378	AS115389
13"	AS115335	AS115346	AS115357	AS115368	AS115379	AS115390
14"	AS115336	AS115347	AS115358	AS115369	AS115380	AS115391
15"	AS115337	AS115348	AS115359	AS115370	AS115381	AS115392
16"	AS115338	AS115349	AS115360	AS115371	AS115382	AS115393
17"	AS115339	AS115350	AS115361	AS115372	AS115383	AS115394
18"	AS115340	AS115351	AS115362	AS115373	AS115384	AS115395
19"	AS115341	AS115352	AS115363	AS115374	AS115385	AS115396
20"	AS115342	AS115353	AS115364	AS115375	AS115386	AS115397



BREATHABLE BACK UPHOLSTERY (CONT.)

TENSION ADJUSTABLE FOR TiLite™ ZR (SERIES 2), ZRA (SERIES 2),
AERO Z (SERIES 2), ZR, ZRA, AERO Z, Z10, TR, YR, BB, AERO R, TR (SERIES 3),
AERO T, TRA (SERIES 2)

TABLE TWO - BACK STRAPS

BACKSTRAP WIDTH	10"-11" CHAIR WIDTH	12"-13" CHAIR WIDTH	14"-15" CHAIR WIDTH	16"-17" CHAIR WIDTH	18"-19" CHAIR WIDTH	20" CHAIR WIDTH
1" WIDE	CP102246	CP101517	CP101518	CP101519	CP101520	CP101521
2" WIDE	CP102232	CP100492	CP100494	CP100496	CP100498	CP101522

*NOTE: FOR ITEM 4

8.5"-11.5" (10"-12.5" FOR PUSH) BACK HEIGHTS USE ONE 1" BACK STRAP AND ONE 2" BACK STRAP

12"-16.5" (12"-16.5" FOR PUSH) BACK HEIGHTS USE TWO 2" BACK STRAPS

16"-20.5" (16"-20.5" FOR PUSH) BACK HEIGHTS USE THREE 2" BACK STRAPS

TABLE THREE - BACKREST EARS

CHAIR WIDTH	STANDARD ADJ. HT 8.5"-11.5"	STANDARD ADJ. HT 12"-16.5"	STANDARD ADJ. HT 16"-20.5"	ADJ. HT 10"-12.5" W/PUSH HANDLES	ADJ. HT 12"-16.5" W/PUSH HANDLES	ADJ. HT 16"-20.5" W/PUSH HANDLES
10"-11" RIGHT	CP102233	CP102234	CP102234	CP102235	CP102236	CP102236
10"-11" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322
12"-13" RIGHT	CP101502	CP100049	CP100049	CP101496	CP100323	CP100323
12"-13" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322
14"-15" RIGHT	CP101504	CP100264	CP100264	CP101497	CP100324	CP100324
14"-15" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322
16"-17" RIGHT	CP101505	CP100265	CP100265	CP101498	CP100325	CP100325
16"-17" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322
18"-19" RIGHT	CP101506	CP100266	CP100266	CP101499	CP100326	CP100326
18"-19" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322
20" RIGHT	CP101507	CP101508	CP101508	CP101500	CP101494	CP101494
20" LEFT	CP101503	CP100050	CP100050	CP101495	CP100322	CP100322

TABLE FOUR - BACK PAD TENSION ADJUSTABLE

CHAIR WIDTH	8.5"-11.5" BACK HEIGHT	12"-16.5" BACK HEIGHT	16"-20.5" BACK HEIGHT
	10"-12.5" BACK HEIGHT W/ PUSH	12"-16.5" BACK HEIGHT W/ PUSH	16"-20.5" BACK HEIGHT W/ PUSH
10"	CP106976	CP106987	CP106998
11"	CP106977	CP106988	CP106999
12"	CP106978	CP106989	CP107000
13"	CP106979	CP106990	CP107001
14"	CP106980	CP106991	CP107002
15"	CP106981	CP106992	CP107003
16"	CP106982	CP106993	CP107004
17"	CP106983	CP106994	CP107005
18"	CP106984	CP106995	CP107006
19"	CP106985	CP106996	CP107007
20"	CP106986	CP106997	CP107008

